

Mast Uri® Well

Order code: SCANURIWELL

Sample preparation and dispensing are critical aspects of the urine screening process.

The Mast Uri® Well is a dispensing aid which minimises sample preparation error. Individual green LED lights, in a 96-well template, automatically indicate the correct well for each sample.



This process is repeated until all the samples have been transferred to the 96-well plate, which becomes the inoculum plate for subsequent inoculation of the Mast Uri® Plate Set(s).

An inoculum plate containing 96 samples may be confidently prepared within 15 minutes.

Specifications	
Weight	0.61Kg
Dimensions	(W x D x H) 100mm x 270mm x 200mm
Power Requirements	5V, 500mA via USB connection
Connectivity	USB



A scan of the specimen barcode will illuminate one of the 96 LED's, indicating the appropriate well for the sample to be dispensed.

