

DiscMaster™

Antibiotic disc dispenser
New and improved



- Smooth ergonomic action
- Compatible with new MastDiscs™ format
- Innovative space-saving storage canister

DiscMaster™ 3

The DiscMaster™ 3 third-generation dispenser provides a reliable and robust solution to dispensing antibiotic susceptibility testing (AST) discs in cartridges within the MastDiscs™ range. The dispenser is designed to complement Mast's new state-of-the-art cartridge manufacturing facilities.

The DiscMaster™ 3 disc dispenser is available in a convenient 6-place format, with an easy-to-use single-handed operation, permitting smooth and rapid integration into routine laboratory workflows.

Ergonomic design

Re-engineered to maintain optimal performance for high-throughput laboratories.

Positive disc placement

AST discs positively placed onto agar surface permitting accurate, centred positioning on each dispensing cycle.

Lock-off indicates empty cartridge

Automatic lock-off alerting user to an empty cartridge.

Height adjustable

Fully height adjustable skirt, accommodating variable depths of culture media.

Audible location of cartridge position

Audible 'click' providing reassurance that cartridges are correctly seated in the dispenser.

Seal-tight canister

Enhanced canister improving rigidity and protection against moisture. **NEW** silica gel capsules with improved colour reaction.

Ordering Information

Order Code	Product	Pack Size
MDD63	DiscMaster™ 3 Dispenser - 6 place	1
SHD1	Single Cartridge Hand Dispenser	1
SHD5	Single Cartridge Hand Dispenser	5
SILICA63	Silica Gel Capsules for DiscMaster™ 3	4

v 1.0 CEA 04/13

United Kingdom
Mast Group Ltd.
Mast House
Derby Road, Bootle
Merseyside L20 1EA

Tel: + 44 (0) 151 933 7277
Fax: + 44 (0) 151 944 1332
e-mail: sales@mastgrp.com

Germany
Mast Diagnostica GmbH
Feldstrasse 20
DE-23858 Reinfeld

Tel: + 49 (0) 4533 2007 0
Fax: + 49 (0) 4533 2007 68
e-mail: mast@mast-diagnostica.de

France
Mast Diagnostic
115 Rue Jules Barni
CS91106
80011 Amiens CEDEX 1

Tél: + 33 (0) 322 80 80 67
Fax: + 33 (0) 322 80 99 22
e-mail: info@mast-diagnostic.fr